X-Safe 240 / X-safe 300EC – single use nitril gloves – sizes S-XL

User information sheet

These gloves are CE-certified in accordance with Regulation EU 2016/425 and were tested to the European Standards EN ISO 21420:2020 “General requirements for working gloves” and EN ISO 374-1 :2016+A1:2018 “Chemical protective gloves” and EN ISO 374-5:2016 “protection against micro organisms.

For testing, certification and annual inspection according to Module C2, the following Notified Body is involved:

NB 0598
SGS Fimko Ltd.
Takomotie 8
FI-00380 Helsinki
Finland

Application:

General: Finger dexterity was tested according to EN ISO 21420:2020 and reached level 5 (Classification minimum level 1 and maximum level 5)

Note to Gloves X-Safe:

With a total length of approximately 240 mm, they are made to fit for special purpose where longer cuffs would empediment certain tasks.

Chemical hazards – EN ISO 374-1:2016+A1:2018

This product standard specifies test methods and performance requirements for gloves which offer protection against dangerous chemicals.

This standard distinguishes 3 test methods; Penetration resistance (water and air leakage), Permeation resistance to liquid chemicals and degradation resistance to liquid chemicals.

Penetration test (water and air leakage) :PASS

The gloves are classified in three types according to the resistance to permeation:

Type A	Permeation breakthrough time \geq 30 minutes for at least 6 products from the list below
Type B	Permeation breakthrough time \geq 30 minutes for at least 3 products from the list below
Type C	Permeation breakthrough time \geq 10 minutes for at least 1 product from the list below

The Type B pictogram is accompanied by a 3-digit code. This code refers to 3 chemicals (from the list of 18 defined chemicals), for which a breakthrough time of at least 30 minutes has been recorded.

EN ISO 374-1:2016+A1:2018/Type B

J K T

Code	Chemical	Cas no.	Category
A	Methanol	67-56-1	Primary alcohol
B	Acetone	67-64-1	Ketone
C	Acetonitrile	75-05-8	Nitrile Compound
D	Dichloromethane	75-09-2	Chlorinated paraffin
E	Carbon disulphide	75-15-0	Sulphur containing organic compound
F	Toluen	108-88-3	Aromatic hydrocarbon
G	Diethylamine	109-89-7	Amine
H	Tetrahydrofuran	109-99-9	Heterocyclic and ethereal compound
I	Ethyl acetate	141-78-6	Ester
J	n-Heptan	142-82-5	Saturated hydrocarbon
K	Sodium hydroxide 40%	1310-73-2	Inorganic base
L	Sulphuric acid 96%	7664-93-9	Inorganic mineral acid
M	Nitric acid 65%	697-37-2	Inorganic mineral acid, oxidizing
N	Acetic acid 99%	64-19-7	Organic acid
O	Ammonia 25%	1336-21-6	Organic base
P	Hydrogene peroxide 30%	7722-84-1	Peroxide
S	Hydrogene fluoride 40%	7664-39-3	Inorganic mineral acid, contact
T	Formaldehyde 37%	50-00-0	Aldehyde

Performance level in accordance with EN ISO 374-1

Measured break through time (min)	Permeation performance level	Performance level
>10	1	
>30	2	
>60	3	n-Heptan – code J
>120	4	
>240	5	
>480	6	Sodium Hydroxide – code K Formaldehyde – code T

Note : The information listed above does not reflect the actual duration of protection at the workplace due to other factors influencing the performance such as temperature, abrasion, degradation, etc.

Performance level in accordance with 374-4:2013 - Resistance to degradation by chemicals

Chemical	Degradation mean	Degradation standard deviation
n-Heptan – code J	47,7%	9,1%
Sodium Hydroxide 40% – code K	2,0%	3,8%
Formaldehyde 37% – code L	33,4%	4,5%

Microbiological hazards – EN ISO 374-5:2016

This product standard specifies test methods and performance requirements for gloves which offer protection against hazardous micro-organisms (bacteria, fungi, viruses).

The applicable test standard is ISO 16604:2004 (method B) and the outcome is “PASS” when no bacteria penetrate the test samples.

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<p>Before usage, inspect the gloves for any defect or imperfections. If in doubt do not use the gloves, get a new pair.</p> <p>Penetration resistance was assessed under laboratory conditions and relates only to the tested specimen.</p>
<p>CAUTION: The levels do not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, abrasion, degradation. Do not use glove which shows signs of imperfection or deterioration on opening the package. Discard immediately gloves which have become worn out or damaged when used. Never wear gloves which are wet or dirty inside to avoid dermatitis.</p>
<p>INGREDIENTS: This product is made from 100% Carboxylated Butadiene Acrylonitrile Copolymer and does not contain natural latex which may cause allergic reactions.</p>
<p>CLEANING: N/A as it is a disposable product.</p>
<p>SHELF LIFE: 5 years from manufacturing date.</p>
<p>EU DECLARATION OF CONFORMITY : Can be found at www.lmtex.dk</p>

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