

X-Soft 300 mm

Nitrilhandske

- En blød og slidstærk nitrilhandske med forbedret elasticitet og større komfort
- 100 % nitril, en syntetisk copolymer
- Finger-teksturen giver et fremragende greb
- Ensartet høj kvalitet
- Høj komfort, god pasform og sensitivitet
- Rullet kantafslutning for bedre holdbarhed
- Modstandsdygtig over for flere kemikalier
- 300 mm længde for bedre sikkerhed
- Tykkelse:
Cuff 0.07-0.08 mm
Palm 0.09-0.10 mm
Finger 0.13-0.16 mm



Varenummer	Størrelse	Farve	Pakning
1514	Small	Blå	100 stk. / æske, 10 æsker / kt.
1515	Medium	Blå	100 stk. / æske, 10 æsker / kt.
1516	Large	Blå	100 stk. / æske, 10 æsker / kt.
1517	X-large	Blå	100 stk. / æske, 10 æsker / kt.
1518	XX-large	Blå	100 stk. / æske, 10 æsker / kt.

EN 420 / EN 374 1-2-3 / EN 455 1-2-3

CE 0120
Directive 89/686/EEC

Kemikalie Test X-Soft 300 Nitril Handske

Kemikalie	CAS NO:	Koncentration	Gennembrudstid
Acetic Acid	64-19-7	10%	480 min
Acetonitrile	75-05-8	99,9 %	1 min
Acrylamide	79-06-1	40 %	480 min
Mixture of 50% Acetonitrile and 2.5% Trifluoroacetic acid in Water	Detection af ACN	Mixture prepared with 99% purity of Acetonitrile and Trifluoroacetic acid	2 min
Mixture of 50% Acetonitrile and 2.5% Trifluoroacetic acid in Water	Detection of TFA	Mixture prepared with 99% purity of Acetonitrile and Trifluoroacetic acid	12 min
Ammonium Hydroxide	1336-21-6	25%	12 min
Boric acid saturated	10043-35-3	98%	480 min
1,4 Butanediol	110-63-4	980%	480 min
Chloroform	67-66-3	98%	0 min
Citric acid	77-92-9	98%	480 min
Cyclohexane	110-82-7	99.9%	83 min
*Dichloromethane	75-09-2	99.0%	1 min
Diethylphthalate	84-66-1	98.0 %	2 min
*Dimethyl Sulfoxide	67-68-5	98.0 %	17 min
Ethanol	64-17-5	99.0%	17 min
*Ethanol	64-17-5	99.0%	14 min
Mixture of Ethanol and Chloroacetic acid	64-17-5 79-11-8	95%	3 min
Ethidium Bromide	1239-45-8	4 %	480 min
Ethylene Glycol	107-21-1	99%	480 min
Ferric Chloride	7705-08-0	99%	480 min
Formaldehyde	50-00-0	37%	480 min
Formamide	75-12-7	98%	480 min
Formic acid	64-18-6	10%	290 min
Hexane	110-54-3	99%	480 min
*Hydrofloric Acid	7664-39-3	48%	3 min
Hydrochloric Acid	7647-01-0	20%	480 min
Hydrogen Peroxide	7722-84-1	20%	280 min
Iso Pentane	78-78-4	99%	75 min

Kemikalie Test X-Soft 300 Nitril Handske

Isopropanol	67-63-0	99.8%	61 min
Methanol	75-05-8	98 %	5 min
Methyl Methacrylate	80-62-6	98 %	1 min
n-Heptane	142-82-5	99%	118 min
Nitric Acid	7697-37-2	10%	480 min
Pentane	109-66-0	98%	42 min
Phenol	108-95-2	85%	8 min
Potassium dichromate	7778-50-9	10%	480 min
* Propanol	71-23-8	99%	10 min
Pyridine	110-86-1	98%	1 min
Sodium Hydroxide	1310-73-2	40 %	480 min
Sodium Hypochlorite, Free Chlorine	7681-52-9	14%	480 min
Sulfuric acid	7664-93-9	98 %	7 min
* Trifluoroacetic Acid	76-05-01	98 %	1 min
Xylene	1330-20-7	98,5 %	3 min

Opdateret 24.10.2018

EN ISO 374-1: 2016 / Type B

EN 16523-1 EN 374-2 EN374-4 EN374-5: 2016

* Dette kemikalie er testet ved temperatur på 33°